

# Clearance Centre

Factory seconds, ex-demo & superseded items

## Elga PURELAB Option R 7 Water Purification Unit™

The Elga PURELAB Option R 7 produces 7L/hr of Type II Pure Water. This is the only recirculating Type II pure water system on the market, ensuring you always have consistent, high specification water for your application.

### Features:

- Optimized resin mixes to increase the ion exchange capacity of consumables and minimize running costs
- Quick and easy sanitization and replacement of consumables to reduce maintenance time
- Designed to be easy to access whether wall or bench mounted with a convenient dispense tap
- The Option R 7 can be used with the wrap-around reservoirs to save space whilst optimizing purity



Product	Condition	Cat. No.	Clearance Price (ex GST)	Available Stock
<b>Product Name</b>				
PURELAB Option R 7 Water Purification Unit		ELGOR007XXM1	\$4590.00	1

### In Australia:

For customer service, call 1300-735-292

To fax an order, use 1800-067-639

To email an order, [ordersau@thermofisher.com](mailto:ordersau@thermofisher.com)

### In New Zealand:

For customer service, call 0800-933-966

To fax an order, use 0800-329-246

To email an order, [ordersnz@thermofisher.com](mailto:ordersnz@thermofisher.com)

**ThermoFisher**  
SCIENTIFIC