

Clearance Centre

Factory seconds, ex-demo & superseded items

Elga PURELAB Ultra Scientific Ultra-Pure Polishing System™

The Elga PURELAB Ultra is the intelligent Type I ultrapure water purification system for your most critical applications. It incorporates many unique and innovative features that set it apart from other systems and guarantees high water purity of 18.2 MΩ-cm.

Features:

- Twin purification pack “PureSure” system guarantees water purity. It stops organics leaching in to the final water purity
- Real-time TOC monitoring ensures confidence in organic quality (Genetic and Analytic models only)
- PureSure system ensures that work is not interrupted by cartridge replacement
- Ongoing bacterial performance guaranteed through complete sanitization of all wetted parts



Product	Condition	Cat. No.	Clearance Price (ex GST)	Available Stock
Product Name				
PURELAB Ultra Scientific Ultra-Pure Polishing System		ELGULXXXSCM2	\$6390.00	1

In Australia:

For customer service, call 1300-735-292

To fax an order, use 1800-067-639

To email an order, ordersau@thermofisher.com

In New Zealand:

For customer service, call 0800-933-966

To fax an order, use 0800-329-246

To email an order, ordersnz@thermofisher.com

ThermoFisher
SCIENTIFIC